



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex pipette tip Qualitix® 10µL, thin graduated tip

€116.48
plus VAT &
Shipping

Product Images



Description

Socorex Qualitix® pipette tips meet the highest quality standards, offer high flexibility and reliable compatibility. The pipette tips are suitable not only for Socorex® pipettes but also for micropipettes of all brands. Different packaging units (bags, rack and the tipfill™ rack refill system), together with a wide range of tips with and without filters, offer enormous versatility.

Product advantages

- Extensive range, high material quality, autoclavable
- Bags, racks and tipfill™ rack refill system
- Guaranteed consistency, batch after batch
- Certified nuclease-, endotoxin- and ATP-free
- Reduced packaging requires less storage space in the lab
- Versatile compatibility with micropipettes of all brands
- autoclavable

Low-binding pipette tips

Low protein binding tips are made of silicone-free resin. Thanks to greatly reduced surface adhesion, they offer maximum liquid recovery. Available in 10, 200 and 1000 µl volume ranges and as sterile or non-sterile racks. Protease-free certified, they are ideally suited for proteomic research and precise handling of viscous samples.

Additional Information

| | |
|---------------------------------------|-------------------------------------|
| No. | SR-0421 |
| Typ Pipettenspitze | Pipettenspitzen low retention |
| Manufacturer (Brand) | Socorex |
| old neoLab article no. | 7-7886 |
| EAN | 4058072369347 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes |
| Volume MAX | 20 µl |
| Volume MIN | 10 µl |
| sterile | Yes |
| Pipette tips properties | with box |
| Packaging type pipette tips (general) | Multiple racks/ trays each with box |
| Packing units tips | 12 x 96 |

