



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Socorex pipette tip Qualitix® 10µL, thin tip

**€133.22**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

**Socorex Qualitix® pipette tips** meet the highest quality standards, offer high flexibility and reliable compatibility. The pipette tips are not only suitable for Socorex® pipettes but also for micropipettes of all brands. Different packaging units (bags, rack and the tipfill™ rack refill system), together with a wide range of tips with and without filters, offer enormous versatility.

Product advantages

- Extensive range, high material quality, autoclavable
- Bags, racks and tipfill™ rack refill system
- Guaranteed consistency, batch after batch
- Certified nuclease-, endotoxin- and ATP-free
- Reduced packaging requires less storage space in the lab
- Versatile compatibility with micropipettes of all brands
- autoclavable

### **Low-binding pipette tips**

Low protein binding tips are made of silicone-free resin. Thanks to greatly reduced surface adhesion, they offer maximum liquid recovery. Available in 10, 200 and 1000 µl volume ranges and as sterile or non-sterile racks. Protease-free certified, they are ideally suited for proteomic research and precise handling of viscous samples.

## Additional Information

---

|                                       |                                     |
|---------------------------------------|-------------------------------------|
| No.                                   | SR-0426                             |
| Typ Pipettenspitze                    | Pipettenspitzen low retention       |
| Manufacturer (Brand)                  | Socorex                             |
| old neoLab article no.                | 7-7894                              |
| EAN                                   | 4058072369392                       |
| Transport temperature                 | Room temperature                    |
| autoclavable at 121°C                 | Yes                                 |
| Volume MAX                            | 20 µl                               |
| Volume MIN                            | 10 µl                               |
| sterile                               | Yes                                 |
| Pipette tips properties               | with box                            |
| Packaging type pipette tips (general) | Multiple racks/ trays each with box |
| Packing units tips                    | 10 x 96                             |

