



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex pipette tip Qualitix® 200µL, oblique graduated tip

€79.87
plus VAT &
Shipping

Product Images



Description

Socorex Qualitix® pipette tips meet the highest quality standards, offer high flexibility and reliable compatibility. The pipette tips are not only suitable for Socorex® pipettes but also for micropipettes of all brands. Different packaging units (bags, rack and the tipfill™ rack refill system), together with a wide range of tips with and without filters, offer enormous versatility.

Product advantages

- Extensive range, high material quality, autoclavable
- Bags, racks and tipfill™ rack refill system
- Guaranteed consistency, batch after batch
- Certified nuclease-, endotoxin- and ATP-free
- Reduced packaging requires less storage space in the lab
- Versatile compatibility with micropipettes of all brands
- autoclavable

Low-binding pipette tips

Low protein binding tips are made of silicone-free resin. Thanks to greatly reduced surface adhesion, they offer maximum liquid recovery. Available in 10, 200 and 1000 µl volume ranges and as sterile or non-sterile racks. Protease-free certified, they are ideally suited for proteomic research and precise handling of viscous samples.

Additional Information

No.	SR-0435
Typ Pipettenspitze	Pipettenspitzen low retention
Manufacturer (Brand)	Socorex
old neoLab article no.	7-8411
EAN	4058072369484
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Volume MAX	200 µl
Volume MIN	200 µl
sterile	No
Pipette tips properties	with box
Packaging type pipette tips (general)	Multiple racks/ trays each with box
Packing units tips	10 x 96

