



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Socorex pipette tip Qualitix® 1000µL, thin tip

**€90.90**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

**Socorex Qualitix® pipette tips** meet the highest quality standards, offer high flexibility and reliable compatibility. The pipette tips are not only suitable for Socorex® pipettes but also for micropipettes of all brands. Different packaging units (bags, rack and the tipfill™ rack refill system), together with a wide range of tips with and without filters, offer enormous versatility.

Product advantages

- Extensive range, high material quality, autoclavable
- Bags, racks and tipfill™ rack refill system
- Guaranteed consistency, batch after batch
- Certified nuclease-, endotoxin- and ATP-free
- Reduced packaging requires less storage space in the lab
- Versatile compatibility with micropipettes of all brands
- autoclavable

### **Low-binding pipette tips**

Low protein binding tips are made of silicone-free resin. Thanks to greatly reduced surface adhesion, they offer maximum liquid recovery. Available in 10, 200 and 1000 µl volume ranges and as sterile or non-sterile racks. Protease-free certified, they are ideally suited for proteomic research and precise handling of viscous samples.

## Additional Information

---

No.	SR-0453
Typ Pipettenspitze	Pipettenspitzen low retention
Manufacturer (Brand)	Socorex
old neoLab article no.	7-8445
EAN	4058072369668
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Volume MAX	1000 µl
Volume MIN	1000 µl
sterile	No
Pipette tips properties	with box
Packaging type pipette tips (general)	Multiple racks/ trays each with box
Packing units tips	6 x 96

