



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex pipette tip Qualitix® refill system Tipfill, colorless tip 10 µL, thin graduated tip

€87.78
plus VAT &
Shipping

Product Images



Description

Socorex Qualitix® pipette tips meet the highest quality standards, offer high flexibility and reliable compatibility. The pipette tips are suitable not only for Socorex® pipettes but also for micropipettes of all brands. Different packaging units (bags, rack and the tipfill™ rack refill system), together with a large selection of tips with and without filters, offer enormous versatility.

Product advantages

- Extensive range, high material quality, autoclavable
- Bags, racks and tipfill™ rack refill system
- Guaranteed consistency, batch after batch
- Certified nuclease-, endotoxin- and ATP-free
- Reduced packaging requires less storage space in the lab
- Versatile compatibility with micropipettes of all brands
- autoclavable

Qualitix® - tipfill™ refill system

Designed for easy rack loading. Filled trays (96 tips) are immediately transferred to a rack without risk of contamination.

The resealable packaging of the tipfill™ refill system reduces storage space, waste volume and shipping costs.

Additional Information

No.	SR-0416
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7880
EAN	4058072369293
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Volume MAX	10 µl
Volume MIN	10 µl
sterile	No
Pipette tips properties	with certification with filter with box RNase-free DNase-free DNA-free ATP-free
Fluid behavior	hydrophobic
Model pipette tip	Socorex filter tips
Type certificate	Batch certificate from the manufacturer
Packaging type pipette tips (general)	Trays for Socorex tipfill system
Packing units tips	13 x 96
for model pipette	Socorex Acura electro Socorex Acura manual

