



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Socorex pipette tip Qualitix® Blue tip, 1000 µL,
autoclavable, slanted tip**

€78.83
plus VAT &
Shipping

Product Images



Description

Socorex Qualitix® pipette tips meet the highest quality standards, offer high flexibility and reliable compatibility. The pipette tips are suitable not only for Socorex® pipettes but also for micropipettes of all brands. Different packaging units (bags, rack and the tipfill™ rack refill system), together with a large selection of tips with and without filters, offer enormous versatility.

Product advantages

- Extensive range, high material quality, autoclavable
- Bags, racks and tipfill™ rack refill system
- Guaranteed consistency, batch after batch
- Certified nuclease-, endotoxin- and ATP-free
- Reduced packaging requires less storage space in the lab
- Versatile compatibility with micropipettes of all brands
- autoclavable

Additional Information

No.	SR-0458
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Socorex
old neolab article no.	7-8450
EAN	4058072369712
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (blue)
Volume MAX	1000 µl
Volume MIN	100 µl
sterile	No
Pipette tips properties	pyrogen-free with certification with box RNase-free DNase-free ATP-free
Fluid behavior	hydrophobic
Model pipette tip	Socorex pipette tips
Type certificate	Batch certificate from the manufacturer
Packaging type pipette tips (general)	Multiple racks/ trays each with box
Packing units tips	6 x 96
for model pipette	Socorex Acura manual Socorex Acura electro

