



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Socorex pipette tip Qualitix® Colorless universal tip, 200µL, oblique, graduated

**€72.59**  
**plus VAT &**  
**Shipping**

### Product Images

---

## Description

---

**Socorex Qualitix® pipette tips** meet the highest quality standards, offer high flexibility and reliable compatibility. The pipette tips are suitable not only for Socorex® pipettes but also for micropipettes of all brands. Different packaging units (bags, rack and the tipfill™ rack refill system), together with a large selection of tips with and without filters, offer enormous versatility.

Product advantages

- Extensive range, high material quality, autoclavable
- Bags, racks and tipfill™ rack refill system
- Guaranteed consistency, batch after batch
- Certified nuclease-, endotoxin- and ATP-free
- Reduced packaging requires less storage space in the lab
- Versatile compatibility with micropipettes of all brands
- autoclavable

## Additional Information

No.	SR-0431
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Socorex
old neoLab article no.	7-8406
EAN	4058072369446
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Volume MAX	200 µl
Volume MIN	1 µl
sterile	No
Pipette tips properties	pyrogen-free with certification with box graduated RNase-free DNase-free ATP-free
Fluid behavior	hydrophobic
Type certificate	Batch certificate from the manufacturer
Packaging type pipette tips (general)	Multiple racks/ trays each with box
Packing units tips	10 x 96
for model pipette	Socorex Acura manual Socorex Acura electro
for pipettes / dispensers of the manufacturers	Brand Eppendorf Gilson Thermo Scientific

