



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex pipette tip Qualitix® refill system Tipfill, universal tip 200µL, slanted, graduated

€64.79
plus VAT &
Shipping

Product Images

Description

Socorex Qualitix® pipette tips meet the highest quality standards, offer high flexibility and reliable compatibility. The pipette tips are not only suitable for Socorex® pipettes but also for micropipettes of all brands. Different packaging units (bags, rack and the tipfill™ rack refill system), together with a large selection of tips with and without filters, offer enormous versatility.

Product advantages

- Extensive range, high material quality, autoclavable
- Bags, racks and tipfill™ rack refill system
- Guaranteed consistency, batch after batch
- Certified nuclease-, endotoxin- and ATP-free
- Reduced packaging requires less storage space in the lab
- Versatile compatibility with micropipettes of all brands
- autoclavable

Qualitix® - tipfill™ refill system

Designed for easy rack loading. Filled trays (96 tips) are immediately transferred to a rack without risk of contamination.

The resealable packaging of the tipfill™ refill system reduces storage space, waste volume and shipping costs.

Additional Information

No.	SR-0432
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Socorex
old neoLab article no.	7-8407
EAN	4058072369453
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Volume MAX	200 µl
Volume MIN	1 µl
sterile	No
Pipette tips properties	pyrogen-free with certification graduated RNase-free DNase-free ATP-free
Fluid behavior	hydrophobic
Packaging type pipette tips (general)	Trays for Socorex tipfill system
Packing units tips	10 x 96
for model pipette	Socorex Acura manual Socorex Acura electro
Type certificate	Batch certificate from the manufacturer
for pipettes / dispensers of the manufacturers	Brand Eppendorf Gilson Thermo Scientific

