



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex pipette tip Qualitix® yellow tip 200µL, sterile, oblique

€82.47
plus VAT &
Shipping

Product Images



Description

Socorex Qualitix® pipette tips meet the highest quality standards, offer high flexibility and reliable compatibility. The pipette tips are suitable not only for Socorex® pipettes but also for micropipettes of all brands. Different packaging units (bags, rack and the tipfill™ rack refill system), together with a large selection of tips with and without filters, offer enormous versatility.

Product advantages

- Extensive range, high material quality, autoclavable
- Bags, racks and tipfill™ rack refill system
- Guaranteed consistency, batch after batch
- Certified nuclease-, endotoxin- and ATP-free
- Reduced packaging requires less storage space in the lab
- Versatile compatibility with micropipettes of all brands
- autoclavable

Additional Information

No.	SR-0439
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Socorex
old neoLab article no.	7-8417
EAN	4058072369521
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	yellow
Volume MAX	200 µl
Volume MIN	1 µl
sterile	Yes
Pipette tips properties	pyrogen-free with certification with box RNase-free DNase-free DNA-free ATP-free
Fluid behavior	hydrophobic
Packaging type pipette tips (general)	Multiple racks/ trays each with box
Packing units tips	10 x 96
for model pipette	Socorex Acura manual Socorex Acura electro
Type certificate	Batch certificate from the manufacturer
for pipettes / dispensers of the manufacturers	Brand Eppendorf Gilson Thermo Scientific

