



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex pipette tips Qualitix® Colorless tip, 350 μ L

€65.31
plus VAT &
Shipping

Product Images



Description

Socorex Qualitix® pipette tips meet the highest quality standards, offer high flexibility and reliable compatibility. The pipette tips are suitable not only for Socorex® pipettes but also for micropipettes of all brands. Different packaging units (bags, rack and the tipfill™ rack refill system), together with a wide range of tips with and without filters, offer enormous versatility.

Product advantages

- Extensive range, high material quality, autoclavable
- Bags, racks and tipfill™ rack refill system
- Guaranteed consistency, batch after batch
- Certified nuclease-, endotoxin- and ATP-free
- Reduced packaging requires less storage space in the lab
- Versatile compatibility with micropipettes of all brands
- autoclavable

Additional Information

No.	SR-0444
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Socorex
old neoLab article no.	7-8431
EAN	4058072369576
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Volume MAX	350 µl
Volume MIN	0.1 µl
sterile	No
Pipette tips properties	with certification RNose-free DNose-free
Fluid behavior	hydrophobic
Packaging type pipette tips (general)	Bag
Packing units tips	10 x 100
for model pipette	Socorex Acura manual Socorex Acura electro
Type certificate	Batch certificate from the manufacturer
for pipettes / dispensers of the manufacturers	Brand Eppendorf Gilson Thermo Scientific

