



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Thermo Scientific Finntip 50-micro, 10x384 pcs/pack

**€623.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Standard pipette tips in racks of the highest quality.

The Finntip Flex System consists of user-friendly, soft flexible high-end tips and a rack with space-saving refill kits .

Due to the particularly smooth surface, approximately 50% higher precision can be achieved, especially when pipetting small volumes. Thanks to particularly thin walls, the soft pipette tips can be attached and ejected with little effort.

The racks are stable and non-slip and can be opened easily with one hand.

The refill system saves space and allows quick insertion of a new matrix with a single click.

## Additional Information

---

No.	PO-0275
Manufacturer (Brand)	Thermo Scientific
old neoLab article no.	3-1476
EAN	4058072338343
Transport temperature	Room temperature
Model pipette tip	Thermo Finntip™

