



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific Finntip 5 ml, 5x54 pcs./pack

€85.60
plus VAT &
Shipping

Product Images



Description

Standard pipette tips in racks of the highest quality.

The Finntip Flex System consists of user-friendly, soft flexible high-end tips and a rack with space-saving refill kits .

Due to the particularly smooth surface, approximately 50% higher precision can be achieved, especially when pipetting small volumes. Thanks to particularly thin walls, the soft pipette tips can be attached and ejected with little effort.

The racks are stable and non-slip and can be opened easily with one hand.

The refill system saves space and allows quick insertion of a new matrix with a single click.

Additional Information

No.	PO-0300
Manufacturer (Brand)	Thermo Scientific
old neoLab article no.	3-1484
EAN	4058072338350
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Volume MAX	5000 µl
Volume MIN	1000 µl
sterile	No
Pipette tips properties	with certification mit box
Fluid behavior	hydrophobic
Designation tip (thermal)	5 ml
Packaging type pipette tips (general)	Multiple racks/ trays each with box
Model pipette tip	Thermo Finntip Flex
Type certificate	Declaration of conformity (standardization)
Packaging type pipette tips (thermal)	Rack
for pipettes / dispensers of the manufacturers	Thermo Scientific

