



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific Finntip Flex 300, 10x96 pcs/pack

€88.80
plus VAT &
Shipping

Product Images



Description

Standard pipette tips of the highest quality.

Soft, flexible high-end tips made of PP.

The particularly smooth surface allows approximately 50% higher precision to be achieved, especially when pipetting small volumes. Thanks to particularly thin walls, the soft pipette tips can be attached and ejected with little effort.

10 boxes à 96 tips

Additional Information

No.	PO-0367
Manufacturer (Brand)	Thermo Scientific
old neoLab article no.	3-1478
EAN	4058072338459
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Volume MAX	300 µl
Volume MIN	5 µl
sterile	No
Pipette tips properties	with certification mit box
Fluid behavior	hydrophobic
Designation tip (thermal)	Flex 300
Packaging type pipette tips (general)	Multiple racks/ trays each with box
Packing units tips	10 x 96
Model pipette tip	Thermo Finntip Flex
Type certificate	Declaration of conformity (standardization)
Packaging type pipette tips (thermal)	Rack
for pipettes / dispensers of the manufacturers	Thermo Scientific

