



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific Finntip Flex 1000, 10x96 pcs/pack

€92.50
plus VAT &
Shipping

Product Images



Description

Standard pipette tips of the highest quality.

Soft, flexible high-end tips made of PP.

The particularly smooth surface allows approximately 50% higher precision to be achieved, especially when pipetting small volumes. Thanks to particularly thin walls, the soft pipette tips can be attached and ejected with little effort.

10 boxes à 96 tips

Additional Information

No.	PO-0372
Manufacturer (Brand)	Thermo Scientific
old neoLab article no.	3-1479
EAN	4058072338466
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Volume MAX	1000 µl
Volume MIN	100 µl
sterile	No
Pipette tips properties	with certification mit box
Fluid behavior	hydrophobic
Designation tip (thermal)	Flex 1000
Packaging type pipette tips (general)	Multiple racks/ trays each with box
Packing units tips	10 x 96
Model pipette tip	Thermo Finntip Flex
Type certificate	Declaration of conformity (standardization)
Packaging type pipette tips (thermal)	Rack
for pipettes / dispensers of the manufacturers	Thermo Scientific

