



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf ep Dualfilter T.I.P.S.® SealMax PCR
clean and sterile, 0.1-10 µL, S, 34 mm, dark gray,
(10 racks x 96 tips)

Special Price

€151.50 was

€202.00

**plus VAT &
Shipping**

Product Images



Description

The purple and white filter of the ep Dualfilter T.I.P.S.® SealMax filter tips combines reliable contamination protection of your pipette from the sample liquid with almost 100% protection from aerosols¹⁾ and biomolecules. This dual filter closes immediately upon contact with the overpipetted sample without inhibiting the valuable material. In addition, droplets and splashes are retained by the white layer; the violet layer provides highly effective protection against sample liquid as a second barrier. For more information and certification on specific sample recovery properties, PCR inhibition and aerosol protection, please visit the Eppendorf website in your country.

Additional Information

No.	VB-2645
Article - USP	Reliable protection against liquids Reliable protection against aerosols ¹⁾ and biomolecules Free of PCR inhibitor additives Protects your pipette against aggressive liquids in case of overpipetting Available in PCR clean/Sterile (Sterile)
Typ Pipettenspitze	Filterspitzen
Manufacturer (Brand)	Eppendorf
EAN	4058072426248
Gross weight	1.04
HSNumber	84799070
Storage temperature	Room temperature
Transport temperature	Room temperature
autoclavable at 121°C	No
Net weight	0.93600 kg
Color	dark gray
Volume MAX	10 µl
Volume MIN	0.1 µl
Country of origin	Germany
sterile	Yes
Customs tariff number	84799070
Length	34 mm
Degree of purity	PCR clean
Packaging type Pipette tips (Eppendorf)	Racks

