



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Eppendorf ep Dualfilter T.I.P.S.® SealMax PCR
clean and sterile, 2-200 µL, 55 mm, yellow, 960
tips (10 racks x 96 tips)**

Special Price

€150.00 was

€200.00

**plus VAT &
Shipping**

Product Images



Description

The purple and white filter of the ep Dualfilter T.I.P.S.® SealMax filter tips combines reliable contamination protection of your pipette from the sample liquid with almost 100% protection from aerosols¹⁾ and biomolecules. This dual filter closes immediately upon contact with the overpipetted sample without inhibiting the valuable material. In addition, droplets and splashes are retained by the white layer; the violet layer provides highly effective protection against sample liquid as a second barrier. For more information and certification on specific sample recovery properties, PCR inhibition and aerosol protection, please visit the Eppendorf website in your country.

Additional Information

No.	VB-2576
Article - USP	Reliable protection against liquids Reliable protection against aerosols ¹⁾ and biomolecules Free of PCR inhibitor additives Protects your pipette against aggressive liquids in case of overpipetting Available in PCR clean/Sterile (Sterile)
Typ Pipettenspitze	Filterspitzen
Manufacturer (Brand)	Eppendorf
EAN	4058072425838
Gross weight	1.244
HSNumber	84799070
Storage temperature	Room temperature
Transport temperature	Room temperature
autoclavable at 121°C	No
Net weight	1.12000 kg
Color	yellow
Volume MAX	200 µl
Volume MIN	2 µl
Country of origin	Germany
sterile	Yes
Customs tariff number	84799070
Length	55 mm
Degree of purity	PCR clean
Packaging type Pipette tips (Eppendorf)	Racks

