



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Eppendorf ep Dualfilter T.I.P.S.® SealMax PCR
clean and sterile, 50-1000 µL, 76 mm, blue, 960
tips (10 racks x 96 tips)**

Special Price

€156.75 was

€209.00

**plus VAT &
Shipping**

Product Images



Description

The purple and white filter of the ep Dualfilter T.I.P.S.® SealMax filter tips combines reliable contamination protection of your pipette from the sample liquid with almost 100% protection from aerosols¹⁾ and biomolecules. This dual filter closes immediately upon contact with the overpipetted sample without inhibiting the valuable material. In addition, droplets and splashes are retained by the white layer; the violet layer provides highly effective protection against sample liquid as a second barrier. For more information and certification on specific sample recovery properties, PCR inhibition and aerosol protection, please visit the Eppendorf website in your country.

Additional Information

| | |
|-----------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| No. | VB-2575 |
| Article - USP | Reliable protection against liquids Reliable protection against aerosols ¹⁾ and biomolecules Free of PCR inhibitor additives Protects your pipette against aggressive liquids in case of overpipetting Available in PCR clean/Sterile (Sterile) |
| Typ Pipettenspitze | Filterspitzen |
| Manufacturer (Brand) | Eppendorf |
| EAN | 4058072425821 |
| Gross weight | 1.87 |
| HSNumber | 84799070 |
| Storage temperature | Room temperature |
| Transport temperature | Room temperature |
| autoclavable at 121°C | No |
| Net weight | 1.68300 kg |
| Color | blue |
| Volume MAX | 1000 µl |
| Volume MIN | 50 µl |
| Country of origin | Germany |
| sterile | Yes |
| Customs tariff number | 84799070 |
| Length | 76 mm |
| Degree of purity | PCR clean |
| Packaging type Pipette tips (Eppendorf) | Racks |

