



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf ep Dualfilter T.I.P.S.® PCR clean and sterile, 0.1-10 µL, S, 34 mm, dark gray, 960 tips (10 racks x 96 tips)

Special Price

€146.25 was

€195.00

plus VAT & Shipping

Product Images



Description

ep Dualfilter T.I.P.S. filter tips are filter tips with contamination protection through a two-phase filter. The two filter layers made of flexible, hydrophobic material fit perfectly into the tip cone and retain almost 100% of all aerosols¹ and biomolecules. This filtering effect is achieved by different, clearly defined pore sizes. The layer facing the sample provides protection against drops, splashes and aerosols. The layer facing the tip cone of the pipette serves as a second barrier against contamination and reliably binds biomolecules. A defined air flow rate ensures complete sample recovery and familiar fast pipetting.

Additional Information

No.	VB-2627
Article - USP	Two filter layers provide dual protection from aerosols ¹⁾ and biomolecules to protect your pipette and sample Free of PCR inhibitor additives for reliable results from downstream processes Available in PCR clean/Sterile (sterile un
Typ Pipettenspitze	Filterspitzen
Manufacturer (Brand)	Eppendorf
VGKL number	106978500
EAN	4058072426071
Gross weight	0.86
HSNumber	84799070
Storage temperature	Room temperature
Transport temperature	Room temperature
autoclavable at 121°C	No
Net weight	0.77400 kg
Color	dark gray
Volume MAX	10 µl
Volume MIN	0.1 µl
Country of origin	Germany
sterile	Yes
Customs tariff number	84799070
Length	34 mm
Degree of purity	PCR clean
Packaging type Pipette tips (Eppendorf)	Racks

