



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf epTIPS Box, Epp.Quality, 0.1-5mL, 120 mm, purple, colorless tips, 24Tips

**€43.60**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

epT.I.P.S.® Box 2.0, Eppendorf Quality™, 0.1 – 5 mL, 120 mm, purple, colorless tips, 1 reusable box × 24 tips. Replaces order no. 0030073177. – Contamination-free transfer of trays into the stable work box – System optimized for use on multichannel pipettes – Color-coded trays for easy tip and pipette identification – Tips can be removed from the refill tray with the pipette – Boxes and trays are autoclavable at 121 °C for later use – epT.I.P.S. Box 2.0 and epT.I.P.S. Set available in Eppendorf Quality

## Additional Information

---

No.	VB-2707
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Eppendorf
EAN	4058072671921
Gross weight	0.51
HSNumber	84799070
Storage temperature	Room temperature
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Net weight	0.45900 kg
Material	Polypropylene (PP)
Volume MAX	5 ml
Volume MIN	0.1 ml
Country of origin	Germany
sterile	No
Customs tariff number	84799070
Length	120 mm
Degree of purity	Eppendorf Quality™
Packaging type Pipette tips (Eppendorf)	Box

