



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf epTIPS Box, Epp.Quality, 0.1-20µL, medium gray, colorless tips, 96Tips

€41.50
plus VAT &
Shipping

Product Images



Description

epT.I.P.S.® Box 2.0, Eppendorf Quality™, 0.1 – 20 µl, medium gray, colorless tips, 1 reusable box × 96 tips. Replaces order no. 0030073037. – Contamination-free transfer of trays into the stable work box – System optimized for use on multichannel pipettes – Color-coded trays for easy tip and pipette identification – Tips can be removed from the refill tray with the pipette – Boxes and trays are autoclavable at 121 °C for later use – epT.I.P.S. Box 2.0 and epT.I.P.S. Set available in Eppendorf Quality

Additional Information

No.	VB-2699
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Eppendorf
EAN	4058072671846
Gross weight	0.28
HSNumber	84799070
Storage temperature	Room temperature
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Net weight	0.25200 kg
Material	Polypropylene (PP)
Volume MAX	20 µl
Volume MIN	0.1 µl
Country of origin	Germany
sterile	No
Customs tariff number	84799070
Degree of purity	Eppendorf Quality™
Packaging type Pipette tips (Eppendorf)	Box

