



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf epTIPS Box, Epp.Quality, 2-200µL, yellow, colorless tips, 96Tips

€36.70
plus VAT &
Shipping

Product Images



Description

epT.I.P.S.® Box 2.0, Eppendorf Quality™, 2 – 200 µL, yellow, colorless tips, 1 reusable box × 96 tips. Replaces order no. 0030073070. – Contamination-free transfer of trays into the stable work box – System optimized for use on multichannel pipettes – Color-coded trays for easy tip and pipette identification – Tips can be removed from the refill tray with the pipette – Boxes and trays are autoclavable at 121 °C for later use – epT.I.P.S. Box 2.0 and epT.I.P.S. Set available in Eppendorf Quality

Additional Information

No.	VB-2701
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Eppendorf
EAN	4058072671860
Gross weight	0.302
HSNumber	84799070
Storage temperature	Room temperature
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Net weight	0.27000 kg
Material	Polypropylene (PP)
Volume MAX	200 µl
Volume MIN	2 µl
Country of origin	Germany
sterile	No
Customs tariff number	84799070
Degree of purity	Eppendorf Quality™
Packaging type Pipette tips (Eppendorf)	Box

