



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf epTIPS Box, Epp.Quality, 50–1,000µL, blue, colorless tips, 96Tips

€36.70
plus VAT &
Shipping

Product Images



Description

epT.I.P.S.® Box 2.0, Eppendorf Quality™, 50 – 1,000 µL, blue, colorless tips, 1 reusable box × 96 tips. Replaces order no. 0030073118. – Contamination-free transfer of trays into the stable work box – System optimized for use on multichannel pipettes – Color-coded trays for easy tip and pipette identification – Tips can be removed from the refill tray with the pipette – Boxes and trays are autoclavable at 121 °C for later use – epT.I.P.S. Box 2.0 and epT.I.P.S. Set available in Eppendorf Quality

Additional Information

No.	VB-2703
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Eppendorf
EAN	4058072671884
Gross weight	0.36
HSNumber	84799070
Storage temperature	Room temperature
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Net weight	0.32400 kg
Material	Polypropylene (PP)
Volume MAX	1000 µl
Volume MIN	50 µl
Country of origin	Germany
sterile	No
Customs tariff number	84799070
Degree of purity	Eppendorf Quality™
Packaging type Pipette tips (Eppendorf)	Box

