



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Eppendorf epT.I.P.S.® Racks, Biopur®, 0.1 – 20 µL,  
40 mm, medium gray, colorless tips, 480 tips (5  
racks × 96 tips).**

**€64.80  
plus VAT &  
Shipping**

## Product Images



## Description

---

With epT.I.P.S. – the original Eppendorf "Totally Integrated Pipetting System" – the tips are optimally matched to the Eppendorf pipettes and, of course, also meet the requirements of EN ISO 8655. This results in minimal attachment and ejection forces with maximum tightness. epT.I.P.S. pipette tips can also be used with pipettes from other manufacturers. Optimal wetting properties, high transparency and special, certified purity grades are the visible expression of our product and production philosophy. epT.I.P.S. Racks Biopur meet the highest demands of the medical, pharmaceutical and food industries. The racks are completely sealed and intended for single use. epT.I.P.S. Racks Biopur are available for volumes from 20 µL to 5 mL.

## Additional Information

No.	VB-2607
Article - USP	Eppendorf Biopur pipette tips in disposable racks have the highest biological purity grade Safe PCR clean, sterile, pyrogen-, ATP- and DNA-free All pre-sterilized epT.I.P.S. and ep Dualfilter T.I.P.S. variants are available in racks<b
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Eppendorf
VGKL number	106975226
EAN	4058072425876
Gross weight	0.546
HSNumber	84799070
Storage temperature	Room temperature
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Net weight	0.49200 kg
diameter	0,8 mm
Color	medium gray
Material	Polypropylene (PP)
Volume MAX	20 µl
Volume MIN	0.1 µl
Country of origin	Germany
sterile	Yes
Customs tariff number	84799070
Length	40 mm
Degree of purity	Biopur
Packaging type Pipette tips (Eppendorf)	Racks

