



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf epTIPS Racks, Biopur®, 20–300µL, 55 mm, orange, colorless tips, 480Tips

€60.30
plus VAT &
Shipping

Product Images



Description

epT.I.P.S.® Racks, Biopur®, 20 - 300 µL, 55 mm, orange, colorless tips, 480 tips (5 racks × 96 tips). Replaces order no. 0030075056 - Eppendorf Biopur pipette tips in disposable racks have the highest biological purity level - Safely PCR clean, sterile, pyrogen-, ATP- and DNA-free - All pre-sterilized epT.I.P.S. and ep Dualfilter T.I.P.S. variants are available in racks - Sterile packaging complies with DIN EN ISO 11607 and DIN EN 868-2-10 standards - Each production batch of racks is continuously tested by an independent laboratory - Batch-specific certificates available at: www.eppendorf.com/certificates - Packaged in racks of 384, 96, 48 or 24 tips - For use in the pharmaceutical and food industries, molecular biology and cell technology - Modern, slim rack with resealable lid - Improved sustainability due to 20% lower polypropylene consumption per rack

Additional Information

No.	VB-2609
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Eppendorf
VGKL number	106975242
EAN	4058072425890
Gross weight	0.65
HSNumber	84799070
Storage temperature	Room temperature
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Net weight	0.58500 kg
diameter	1 mm
Material	Polypropylene (PP)
Volume MAX	300 µl
Volume MIN	20 µl
Country of origin	Germany
sterile	Yes
Customs tariff number	84799070
Length	55 mm
Degree of purity	Biopur
Packaging type Pipette tips (Eppendorf)	Racks

