



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Eppendorf epT.I.P.S.® Reloads, Eppendorf
Quality™, 0.1 – 10 µL, 34 mm, dark gray, colorless
tips, 960 tips (10 trays × 96 tips)**

**€83.20
plus VAT &
Shipping**

Product Images



Description

With epT.I.P.S.®-, the original Eppendorf "Totally Integrated Pipetting System"-, the tips are optimally matched to the Eppendorf pipettes and, of course, also meet the requirements of EN ISO 8655. This results in minimal attachment and ejection forces with maximum tightness. epT.I.P.S. pipette tips can also be used with pipettes from other manufacturers. Optimal wetting properties, high transparency and special, certified purity grades are the visible expression of our product and production philosophy. The epT.I.P.S. Reloads in the Eppendorf Quality are high-quality original Eppendorf pipette tips in refill trays and part of the refill system - box and reload. The refill trays are autoclavable. They are available for volumes from 10 µL to 2.5 mL.

Additional Information

No.	VB-0157
Article - USP	Less waste compared to disposable racks The refill trays are autoclavable at up to 121 °C in the packaging and without film for later use Depending on the tip size, the trays in the refill system are placed against each other or over each other
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Eppendorf
VGKL number	106973371
EAN	4058072404956
Gross weight	0.5
HSNumber	84799070
Storage temperature	Room temperature
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Net weight	0.45000 kg
Color	dark gray
Material	Polypropylene (PP)
Volume MAX	10 µl
Volume MIN	0.1 µl
Country of origin	Germany
sterile	No
Customs tariff number	84799070
Length	34 mm
Degree of purity	Eppendorf Quality™
Packaging type Pipette tips (Eppendorf)	Reloads

