



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Eppendorf epT.I.P.S.® Reloads, Eppendorf
Quality™, 20 – 300 µL, 55 mm, orange, colorless
tips, 960 tips (10 trays × 96 tips)**

**€71.20
plus VAT &
Shipping**

Product Images



Description

With epT.I.P.S.®, the original Eppendorf "Totally Integrated Pipetting System"-, the tips are optimally matched to the Eppendorf pipettes and, of course, also meet the requirements of EN ISO 8655. This results in minimal attachment and ejection forces with maximum tightness. epT.I.P.S. pipette tips can also be used with pipettes from other manufacturers. Optimal wetting properties, high transparency and special, certified purity grades are the visible expression of our product and production philosophy. The epT.I.P.S. Reloads in the Eppendorf Quality are high-quality original Eppendorf pipette tips in refill trays and part of the refill system - box and reload. The refill trays are autoclavable. They are available for volumes from 10 µL to 2.5 mL.

Additional Information

No.	VB-0165
Typ Pipettenspitze	Pipettenspitzen
Article - USP	Less waste compared to disposable racks The refill trays are autoclavable at up to 121 °C in the packaging and without film for later use Depending on the tip size, the trays in the refill system are placed against each other or over each other
Manufacturer (Brand)	Eppendorf
VGKL number	106973452
EAN	4058072404994
Gross weight	0.9
Transport temperature	Room temperature
HSNumber	84799070
autoclavable at 121°C	Yes
Storage temperature	Room temperature
Net weight	0.81000 kg
Color	orange
Material	Polypropylene (PP)
Volume MAX	300 µl
Volume MIN	20 µl
sterile	No
Country of origin	Germany
Customs tariff number	84799070
Length	55 mm
Degree of purity	Eppendorf Quality™
Packaging type Pipette tips (Eppendorf)	Reloads

