



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Eppendorf epT.I.P.S.® Reloads, PCR clean, 50 -
1,250 µL L, 103 mm, dark green, colorless tips,
960 tips (10 trays × 96 tips)**

**€97.70
plus VAT &
Shipping**

Product Images



Description

With epT.I.P.S.®, the original Eppendorf "Totally Integrated Pipetting System"-, the tips are optimally matched to the Eppendorf pipettes and, of course, also meet the requirements of EN ISO 8655. This results in minimal attachment and ejection forces with maximum tightness. epT.I.P.S. pipette tips can also be used with pipettes from other manufacturers. Optimal wetting properties, high transparency and special, certified purity grades are the visible expression of our product and production philosophy. The epT.I.P.S. Reloads PCR clean are high-quality original Eppendorf pipette tips in refill trays and part of the refill system - box and reload. The refill trays are autoclavable. They are available for volumes from 10 µL to 2.5 mL.

Additional Information

| | |
|-----------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| No. | VB-0176 |
| Article - USP | Less waste compared to disposable racks The refill trays are autoclavable at up to 121 °C in the packaging and without film for later use Depending on the tip size, the trays in the refill system are placed against each other or over each other |
| Typ Pipettenspitze | Pipettenspitzen |
| Manufacturer (Brand) | Eppendorf |
| VGKL number | 106973649 |
| EAN | 4058072405069 |
| Gross weight | 1.55 |
| HSNumber | 84799070 |
| Storage temperature | Room temperature |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes |
| Net weight | 1.39500 kg |
| Color | dark green |
| Material | Polypropylene (PP) |
| Volume MAX | 1250 µl |
| Volume MIN | 50 µl |
| Country of origin | Germany |
| sterile | No |
| Customs tariff number | 84799070 |
| Length | 103 mm |
| Degree of purity | PCR clean |
| Packaging type Pipette tips (Eppendorf) | Reloads |

