



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Eppendorf epT.I.P.S.® Singles, Biopur®, 50 – 1,000  
µL, 71 mm, blue, colorless tips, 100 tips,  
individually packed**

**€33.60  
plus VAT &  
Shipping**

## Product Images



## Description

---

With epT.I.P.S.®-, the original Eppendorf "Totally Integrated Pipetting System"-, the tips are optimally matched to the Eppendorf pipettes and, of course, also meet the requirements of EN ISO 8655. This results in minimal attachment and ejection forces with maximum tightness. epT.I.P.S. pipette tips can also be used with pipettes from other manufacturers. Optimal wetting properties, high transparency and special, certified purity grades are the visible expression of our product and production philosophy. epT.I.P.S. Singles are individually packaged tips in the most comprehensive purity grade Eppendorf Biopur. These tips are available for volumes of 20 µL, 200 µL and 1,000 µL.

## Additional Information

No.	VB-0096
Article - USP	Eppendorf Biopur pipette tips, individually packed Safe PCR clean, sterile, pyrogen-, ATP- and DNA-free For use in pharmaceutical and food industry, molecular biology and cell technology Lot number and expiry date on each B
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Eppendorf
EAN	4058072404499
Gross weight	0.205
HSNumber	84799070
Storage temperature	Room temperature
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Net weight	0.18500 kg
Color	blue
Material	Polypropylene (PP)
Volume MAX	1000 µl
Volume MIN	50 µl
Country of origin	Germany
sterile	Yes
Customs tariff number	84799070
Length	71 mm
Degree of purity	Biopur
Packaging type Pipette tips (Eppendorf)	Singles

