



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® neoMultiTray pipette tips 200-1000 µl natural, 16 x 60 pcs./pack

**€26.53**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

In the neoMultiTray plus, the stock pack, the pipette tips are placed on plastic pallets. For clear volume identification, the neoMultiTrays are provided with colored frames.

Each pallet can be removed individually and securely snapped into the neoMultibox plus (Order No. 7-3040) with a light click.

The German-made pipette tips are made of high-quality PP, have a tight fit thanks to a precisely fitting cone, excellent transparency and outstanding volume accuracy.

They fit most standard microliter pipettes, e.g. Eppendorf, Gilson, Socorex, Finn, Biohit, neoPette.

Accessories: neoMultibox plus (Order No. 7-3040)

## Additional Information

---

No.	7-3043
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	neoLab
EAN	4058072038557
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Volume MAX	1000 µl
Volume MIN	200 µl
sterile	No
Pipette tips properties	with certification mit box
Fluid behavior	hydrophobic
Model pipette tip	neoLab neoMultiTray Tips
Type certificate	Declaration of conformity (standardization)
Packaging type pipette tips (general)	Refill. stackable
Packing units tips	16 x 60
for pipettes / dispensers of the manufacturers	Eppendorf Gilson Brand Thermo Scientific neoLab

