



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Pipetting Aids, Easypet 3 incl. power supply, wall holder, storage stand, 2 membrane filters 0.45µm (non-sterile)

**€544.00
plus VAT &
Shipping**

Product Images



Description

Eppendorf Easypet® 3, incl. power supply, wall holder, storage stand, 2 membrane filters 0.45 µm, 0.1 – 100 mL

Experience a new dimension in electronic pipetting with the Eppendorf Easypet® 3. State-of-the-art technological developments give you full speed control and highest precision at the same time. When developing the lightweight, ergonomic Easypet 3, particular attention was paid to ensuring that it fits comfortably in the hand and adapts to its shape for fatigue-free pipetting. Intuitive and convenient speed adjustment simply at the touch of a finger. Lightweight, ergonomic design for fatigue-free pipetting. Easy-to-read, illuminated LED display of remaining battery life for high reliability. Lithium polymer battery for long, cable-free operation. Continuous pump speed adjustment. Continued operation possible while charging. Autoclavable pipette adapter for sterile applications. Pick-up cone with quick coupling facilitates replacement of membrane filters.

Additional Information

| | |
|--------------------------|--|
| No. | VB-1007 |
| Manufacturer (Brand) | Eppendorf |
| EAN | 4058072171469 |
| Transport temperature | Room temperature |
| Color | gray |
| suitable for | Graduated pipettes |
| Pipette Helper Features | with charging station membrane filter 0.45µm |
| Energy supply | Power supply unit |
| Functionality | electronic |
| for pipettes with volume | up to 100 ml |
| Type pipetting aid | electric |

