



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Brand Micro pipette controller, for pipettes up to 1 ml

€49.80
plus VAT &
Shipping

Product Images



Description

The pipette controller micro is an indispensable accessory for sampling with disposable micropipettes with ring mark and many other small-volume pipettes up to 1 ml (e.g. blood mixing and blood glucose pipettes). Steam sterilization resistant (121 °C). The micro pipetting assistant reduces the risk of infection. Contaminated disposable micropipettes up to 50 µl can be discarded using the ejector integrated in the device. This largely prevents the transmission of dangerous viruses such as hepatitis B or HIV. The micro pipetting controller is extremely light and fits comfortably in the hand. Because of its ergonomic shape and easy handling, the micro-classic pipetting controller offers the comfort needed for this strenuous work. Especially in the IVF laboratory, but also in the medical laboratory, it proves to be a reliable partner. Suitable for disposable micropipettes with ring mark and other small volume pipettes up to 1 ml (e.g. blood mixing pipettes). The adapter can be rotated to the left and right. Adapter and suction tubing are steam sterilization resistant at 121 °C. The micro-classic reduces the risk of infection when handling infectious material.

Additional Information

No.	BN-0662
Manufacturer (Brand)	Brand
VGKL number	274925800
EAN	4058072183417
Gross weight	0.1236
HSNumber	84798997
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.10560 kg
Country of origin	Germany
Packing width	0.07400 m
Packing height	0.03100 m
Packing depth	0.20200 m
Packaging volume	463.38800000 cc
Customs tariff number	84798997

