



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Thermo Scientific S1 pipetting aid

**€560.00**  
**plus VAT &**  
**Shipping**

### Product Images

---

## Description

---

S1 Pipetting Aid The Thermo Scientific™ S1 Pipetting Aid is a lightweight, cordless pipetting aid that provides longer, fatigue-free pipetting—for use with glass or plastic graduated and volumetric serological pipettes. Current rating: 0.2A Battery type: lithium-ion Electrical requirements: 100 to 240 V, 50/60 Hz Filter type: 0.45 micron Standard For use with (devices): All major brands of glass or plastic serological pipettes, 1 to 100 ml Frequency: 50/60 Hz Includes: S1 pipetting aid, universal charger, table stand, holder for wall mounting, 1 ml pipette tray, user manual and warranty card Charging time: 3 hours Voltage: 100 to 240 V Warranty: 2 years if registered on our website Dispensing speed: 8 (and dispensing by gravity) For use with: Graduated and volumetric glass serological pipette, Graduated and volumetric plastic serological pipette Volume (metric): 1 to 100 ml Unit Size: Each

## Additional Information

---

No.	PO-0465
Manufacturer (Brand)	Thermo Scientific
EAN	4066292158368
Gross weight	0.93
Transport temperature	Room temperature
Color	White
Country of origin	China
Packaging volume	0.00600000 ccm
Customs tariff number	84798997

