



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific S1 pipetting aid

€560.00
plus VAT &
Shipping

Product Images

Description

S1 Pipetting Aid The Thermo Scientific™ S1 Pipetting Aid is a lightweight, cordless pipetting aid that provides longer, fatigue-free pipetting—for use with glass or plastic graduated and volumetric serological pipettes. Current: 0.2 A Battery type: lithium-ion Electrical requirements: 100 to 240 V 50/60 Hz Filter type: 0.45 micron, standard For use with (devices): Glass or plastic serological pipettes of all major brands, 1 to 100 ml Frequency: 50/60 Hz Includes: S1 pipetting aid, universal charger, table stand, holder for wall mounting, 1 ml pipette tray, user manual and warranty card Charging time: 3 hours Voltage: 100 to 240 V Warranty: 2 years with web registration Dispensing speed: 8 (+ dispensing by gravity) For use with: Graduated and volumetric glass serological pipette, Graduated and volumetric plastic serological pipette Volume (metric): 1 to 100 ml Unit Size: Each

Additional Information

No.	PO-0474
Manufacturer (Brand)	Thermo Scientific
EAN	4066292158450
Gross weight	0.928
Transport temperature	Room temperature
Color	blue
Country of origin	China
Packaging volume	0.00600000 ccm
Customs tariff number	84798997

