



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Thermo Scientific S1 pipetting aid

**€560.00**  
**plus VAT &**  
**Shipping**

### Product Images

---

## Description

---

S1 Pipetting Aid The Thermo Scientific™ S1 Pipetting Aid is a lightweight, cordless pipetting aid that provides longer, fatigue-free pipetting—for use with glass or plastic graduated and volumetric serological pipettes. Battery type: lithium-ion Dispensing speed: Unshielded ; Tube: 13 For use with: Graduated and volumetric glass serological pipette, Graduated and volumetric plastic serological pipette Volume (metric): 1 to 100 ml Filter type: 0.45 micrometer, standard For use with (devices): Glass or plastic serological pipettes of all major brands, 1 to 100 ml Includes: S1 pipetting aid, universal charger, bench stand, wall mount bracket, 1 ml pipette tray, user manual and warranty card Recharge time: 3 hours Warranty: 2 years with web registration Electrical requirements: 100 to 240 V 50/60 Hz Voltage: 100 to 240 V Amperage: 0.2 A Frequency: 50/60 Hz Unit Size: Each

## Additional Information

---

No.	PO-0476
Manufacturer (Brand)	Thermo Scientific
EAN	4066292158474
Gross weight	0.93
Transport temperature	Room temperature
Color	green
Country of origin	China
Packaging volume	0.00600000 ccm
Customs tariff number	84798997

