



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Thermo Scientific S1 pipetting aid

**€560.00**  
**plus VAT &**  
**Shipping**

### Product Images

---

## Description

---

S1 Pipetting Aid The Thermo Scientific™ S1 Pipetting Aid is a lightweight, cordless pipetting aid that provides longer, fatigue-free pipetting—for use with glass or plastic graduated and volumetric serological pipettes. Current: 0.2 A Battery type: lithium-ion Electrical requirements: 100 to 240 V 50/60 Hz Filter type: 0.45 micron, standard For use with (devices): Glass or plastic serological pipettes of all major brands, 1 to 100 ml Frequency: 50/60 Hz Includes: S1 pipetting aid, universal charger, table stand, holder for wall mounting, 1 ml pipette tray, user manual and warranty card Charging time: 3 hours Voltage: 100 to 240 V Warranty: 2 years with web registration Dispensing speed: 8 (+ dispensing by gravity) For use with: Graduated and volumetric glass serological pipette, Graduated and volumetric plastic serological pipette Volume (metric): 1 to 100 ml Unit Size: Each

## Additional Information

---

|                       |                     |
|-----------------------|---------------------|
| No.                   | PO-0469             |
| Manufacturer (Brand)  | Thermo Scientific   |
| EAN                   | 4066292158405       |
| Gross weight          | 0.924               |
| Transport temperature | Room temperature    |
| Color                 | transparent (clear) |
| Country of origin     | China               |
| Packaging volume      | 0.00600000 ccm      |
| Customs tariff number | 84798997            |

