



neoLab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neoLab® Line pipette controller, yellow

**€149.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

With the new neoPipettor pipetting aid, you have pipetting safely under control. Comfortable working with glass or plastic pipettes from 0.1 ml to 100 ml volume. Replaceable 0.45 µm membrane filters to minimize the risk of contamination. The speed for liquid aspiration and dispensing can be continuously adjusted.

- Ergonomically shaped handle with smooth-running control knobs
- Comfortable and effortless pipetting
- Low weight only 200 g
- Stepless control knobs for precise and controlled pipetting
- Modern design
- Cordless operation (8 hours without recharging)
- Fast recharge time within 2-3 hours
- Quiet high performance motor
- Stepless speed adjustment
- Soft ribbed silicone adapter for easy pick up of pipettes

Scope of delivery:

- Charger
- Table stand
- Lithium-ion battery

**Product video:**

## Additional Information

No.	DK-0111
Manufacturer (Brand)	neoLab
old neoLab article no.	D-6018
EAN	4058072981228
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	yellow
suitable for	Pipettes
Pipette Helper Features	with infinitely variable speed control with membrane filter with charging station
Number of channels	1
Energy supply	Battery
Spare part	D-6028 D-6026 D-6029 D-6027
Functionality	electric
Type pipetting aid	electric
for pipettes with volume	up to 100 ml

