



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® disposable 8-channel transfer pipette, non-sterile, 25 pcs./pack

€108.20
plus VAT &
Shipping

Product Images



Description

Single-use 8-channel transfer pipette

Quickly and easily wash microplates or transfer media with the unique single-use 8-channel transfer pipette. Without the use of pipette tips or pipettors, the simplest pipetting can be performed.

Sampling or filling of microtiter plates without much effort.

600µl can be pipetted per channel.

Non-sterile.

Additional Information

No.	7-4900
Manufacturer (Brand)	neoLab
EAN	4058072041281
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (crystal clear)
Material	Low-density polyethylene (LDPE)
Volume	0.6 ml
suitable for	Microtiter plates
Number of channels	8
Functionality	mechanical

