




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Brand macro pipetting aid, up to 200 ml, gray

**€62.35**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Mechanical pipette controller for bulb and graduated pipettes from 0.1 to 200 ml.

Simple operation via one-hand function elements.

Low weight of 106 g.

Ideal even for inexperienced users.

A hydrophobic membrane filter protects against liquid penetration.

**Completely autoclavable up to 121°C.**

Accessories: Replacement membrane filter 3 µm (Order No. 1-2601)

Replacement suction bellows (Order No. 1-2602)

Spare valve system (Order No. 1-2603)

## Additional Information

---

No.	BN-0675
Manufacturer (Brand)	Brand
old neolab article no.	1-2600
EAN	4058072167783
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	gray
Material	Polyamide (PA. non-specific)
suitable for	Graduated pipettes full pipettes
Pipette Helper Features	with membrane filter
Number of channels	1
Functionality	mechanical
for pipettes with volume	0.1 to 200 ml
Spare part	BN-0671 BN-0680 BN-0673

