



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Brand macro pipetting aid, up to 200 ml, magenta

€62.35
plus VAT &
Shipping

Product Images



Description

Mechanical pipette controller for bulb and graduated pipettes from 0.1 to 200 ml.

Simple operation via one-hand function elements.

Low weight of 106 g.

Ideal even for inexperienced users.

A hydrophobic membrane filter protects against liquid penetration.

Completely autoclavable up to 121°C.

Accessories: Replacement membrane filter 3 µm (Order No. 1-2601)

Replacement suction bellows (Order No. 1-2602)

Spare valve system (Order No. 1-2603)

Additional Information

No.	BN-0678
Manufacturer (Brand)	Brand
old neolab article no.	1-2608
EAN	4058072167820
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	pink
Material	Polyamide (PA. non-specific)
suitable for	Graduated pipettes full pipettes
Pipette Helper Features	with membrane filter
Number of channels	1
Functionality	mechanical
for pipettes with volume	0.1 to 200 ml
Spare part	BN-0671 BN-0680 BN-0673

