



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Vitlab pipetting aid green, 10 ml

€19.42 As low as
€17.65
**plus VAT &
Shipping**

Product Images



Description

Particularly suitable for pipettes calibrated to the tip.

- aspiration of the medium by turning a knurled wheel
- Quick emptying by pressing a small lever
- soft, elastic flange with threaded collar for secure attachment of the pipettes
- easy to dismantle and clean
- Sizes coded according to colors

Additional Information

No.	LV-0492
Manufacturer (Brand)	Vitlab
old neoLab article no.	6-2301
EAN	4058072292744
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	green
Material	Polypropylene (PP)
suitable for	Graduated pipettes full pipettes
Number of channels	1
Functionality	mechanical
for pipettes with volume	up to 10 ml

