



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Vitlab pipetting aid red, 25 ml

**€20.03** As low

as

**€18.21**

**plus VAT &  
Shipping**

## Product Images

---



## Description

---

Particularly suitable for pipettes calibrated to the tip.

- aspiration of the medium by turning a knurled wheel
- Quick emptying by pressing a small lever
- soft, elastic flange with threaded collar for secure attachment of the pipettes
- easy to dismantle and clean
- Sizes coded according to colors

## Additional Information

---

No.	LV-0493
Manufacturer (Brand)	Vitlab
old neoLab article no.	6-2302
EAN	4058072292751
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	red
Material	Polypropylene (PP)
suitable for	Graduated pipettes full pipettes
Number of channels	1
Functionality	mechanical
for pipettes with volume	up to 25 ml

