



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Howorka pipetting ball 20 ml, autoclavable

**€22.30**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Universally applicable laboratory auxiliary device.

Specially designed opening for immediate sealing effect with all common tube diameters.

No unnecessary valves, disturbing air cushions or rigid connections.

- aspirating and blowing out pipettes, burettes, capillaries
- pipetting and filtering of organic solvents and viscous substances

Easy to clean, autoclavable.

Recommended by the German Employer's Liability Insurance Association as "safe work equipment" .

## Additional Information

---

No.	3-1100
Manufacturer (Brand)	neoLab
EAN	4058072012267
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	Black
Material	Rubber. natural
Volume	20 ml
Number of valves / cocks	0

