



neoLab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neoLab® hand operator set in stand

**€310.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The neoVAQ handle set for aspirators, including a practical stand and various attachments, is the ideal helper for all aspiration problems in cell culture.

Ergonomically shaped handle with pressure-sensitive regulation knob.

Easy and fast exchange of adapters by means of quick coupling on the handle.

### Highlights

- multifunctional use
- Adapters for different applications
- autoclavable, except for the rack
- Material: POM & stainless steel (attachments & handle)
- Stand made of ABS
- Vacuum on/off switch
- connectable to different vacuum sources
- adapted to laboratory vessels & cell culture products (e.g. microtiter plates)

Scope of delivery:

Handle, color light blue

1-channel tip attachment for pipette tips or serological pipettes

1-channel stainless steel attachment 40 mm long

1-channel stainless steel attachment 80 mm long

8-channel stainless steel attachment 40 mm long

8-channel tip attachment with ejector for pipette tips

1-channel tip adapter with ejector

1-channel tip adapter (150 mm)

Stand (ABS), non-autoclavable

## Additional Information

No.	7-0244
Manufacturer (Brand)	neoLab
EAN	4058072158095
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	Blue White
Material	Polyoxymethylene (POM)
suitable for	Extraction systems
Pipette Helper Features	autoclavable
Number of channels	1
Functionality	mechanical
Number of valves / cocks	1
Scope of delivery	Handle, color light blue 1-channel tip attachment for pipette tips or serological pipettes 1-channel stainless steel tip 40 mm long 1-channel stainless steel tip attachment

