



Hirschmann® pipetus®-junior pipetting aid

Product Images



Description

Manual pipetting aid for all pipettes between 0.5 and 25 ml.
Simple and safe, accurate and inexpensive.
One complete stroke aspirates about 2 ml of liquid.
Pressing the OUT button causes the pipette to run out, the intensity of the button pressure determines the run-out speed.
The standard hydrophobic PTFE membrane filter acts as an integrated overfill protection.
Weight: 140 g

Additional Information

No.	NH-1692
Manufacturer (Brand)	Hirschmann
old neoLab article no.	3-9902
EAN	4058072336790

