



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Cryomaster® cabinet insert short combi for 10x136mm and 5x137mm boxes

€233.50
plus VAT &
Shipping

Product Images



Description

Cryomaster cabinet insert short combi
for 10 boxes 136 mm
and for 5 boxes 137 mm

Additional Information

| | |
|----------------------------------|--|
| No. | 4-6214 |
| Manufacturer (Brand) | Cryomaster |
| EAN | 4058072157814 |
| Transport temperature | Room temperature |
| Color | silver |
| Material | Stainless steel |
| Wide | 140 mm |
| Height | 249 mm |
| Length | 700 mm |
| Depth | 700 mm |
| Number of compartments in height | 2 |
| Number of compartments in depth | 5 |
| suitable for cooling in | Ultra Deep Freezer Freezer Refrigerator |
| Type slide-in | Cabinet drawers for 50-boxes Cabinet drawers for 100-boxes |

