



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Cryomaster® cabinet drawer for 72MTP 86x128x18 mm, 6 deep, 12 high

€296.80
plus VAT &
Shipping

Product Images



Description

Cryomaster cabinet drawer for 72 MTP; 86 x 128 x 18 mm
6 deep, 12 high

Additional Information

| | |
|----------------------------------|---|
| No. | 4-6250 |
| Manufacturer (Brand) | Cryomaster |
| EAN | 4058072191757 |
| Transport temperature | Room temperature |
| Color | silver |
| Material | Stainless steel |
| Wide | 138 mm |
| Height | 314 mm |
| Depth | 55 mm |
| Number of compartments in height | 12 |
| Number of compartments in depth | 4 |
| Capacity | 72 |
| Type slide-in | Cabinet inserts for 50 boxes |
| suitable for cooling in | Ultra Deep Freezer Freezer Refrigerator |

