



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Cryomaster® cabinet drawer short for 16 cryoboxes 136x136x103 mm

Special Price

€149.00 was

€233.50

plus VAT &
Shipping

Product Images



Description

Cryomaster® stainless steel cabinet drawers for all boxes with base area 136 x 136 mm (W x D)

Additional Information

No.	4-6210
Manufacturer (Brand)	Cryomaster
EAN	4058072161255
Transport temperature	Room temperature
Color	silver
Material	Stainless steel
Wide	140 mm
Height	414 mm
Length	563 mm
Depth	563 mm
Number of compartments in height	4
Number of compartments in depth	4
Capacity	16
Maximum width box	136 mm
Maximum height box	103 mm
Maximum length box	136 mm
Type slide-in	Cabinet drawers for 100 boxes
suitable for cooling in	Ultra Deep Freezer Freezer Refrigerator

