



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Cryomaster® cabinet drawer for 25 cryoboxes 136x136x53 mm

€233.50
plus VAT &
Shipping

Product Images



Description

Cryomaster® stainless steel cabinet drawers for all boxes with base area 136 x 136 mm (W x D)

Additional Information

No.	4-6212
Manufacturer (Brand)	Cryomaster
EAN	4058072157791
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	silver
Material	Stainless steel
sterile	No
Wide	140 mm
Height	288 mm
Length	563 mm
Depth	700 mm
Number of compartments in height	5
Number of compartments in depth	5
suitable for cooling in	Ultra Deep Freezer Freezer Refrigerator
Maximum width box	136 mm
Maximum height box	53 mm
Maximum length box	136 mm
Type slide-in	Cabinet inserts for 50 boxes

