



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Cryomaster® cabinet drawer w. stop f. 3x4 neoBox-81, 141 x 434 x 245 mm

**Price on
request
plus VAT &
Shipping**

Product Images



Description

For all boxes with outer dimensions 133 x 133 x 50 mm (W x D x H).

Specially manufactured for the neoBox-81 (order no. 2-2910 ff), but also suitable for boxes no. 7-8041 ff and 2-5026

Basic dimension: 141 mm width

Additional Information

No.	2-8493
Manufacturer (Brand)	Cryomaster
EAN	4058072010331
Transport temperature	Room temperature
Color	silver
Material	Stainless steel
Wide	434 mm
Width/inside	410 mm
Height	245 mm
Height inside	55.0 mm
Length	141 mm
Depth	141 mm
suitable for	Boxes/boxes
Length inside	135 mm
Racks properties	with pullout stop
Number of compartments in height	4
Number of compartments in depth	3
Number of fans length	1
suitable for cooling in	Refrigerator Freezer Ultrat freezer Deep freezer Nitrogen container
Maximum width box	133 mm
Maximum height box	50 mm
Maximum length box	133 mm
Type slide-in	Cabinet inserts for neoBox-81

