



Thermo Scientific CryoBox™ freezer racks

Product Images

Description

CryoBox™ Freezer Racks Conveniently store standard CryoBoxes with these freezer racks made of durable stainless steel or aluminum that fit most standard freezers and chest freezers. Construction/Design: Horizontal Dimensions (L x W x H): 22.5 x 44.4 x 14.3 cm Capacity: 12 CryoBoxes/Rack Quantity: 2/cs. Type: Cryogenic Box Storage Rack, Horizontal Unit Size: Case of 2

Additional Information

No.	CD-1539
Manufacturer (Brand)	Thermo Scientific
VGKL number	238750381
EAN	4066292478244
Gross weight	6.76
Country of origin	USA
Packaging volume	0.04800000 ccm
Customs tariff number	73269098

